

**Titan III Wireless Cellular Gateway 4G Introduction
For US models only**

Dear Customer,

Thank you for choosing TITAN III Wireless Cellular Gateway 4G. It has been designed to provide you easy to use and reliable wireless connectivity for data communications.

Before your Gateway can be configured for use with the LIFENET® System, you must associate an active cellular data plan with each Gateway.

If you have purchased the data plan from Stryker/Physio-Control, your Gateway will be shipped activated and your Stryker/Physio-Control Implementation Project Manager will contact you to complete the configuration and registration of your Gateway with the LIFENET System.

If you will be using your own data plan, please follow the instructions below to complete the activation of your Gateway on the cellular network.

The following information will be required to activate the Gateway with your cellular Carrier:

For activation on the Verizon LTE network:

- 1. Find the International Mobile Equipment Identity (IMEI) for your cellular Gateway. For your convenience, this can be found on the peelable portion of the box label, as well as on the label applied to the device. This IMEI # is made up of 15 decimal digits.**
- 2. Contact your local Carrier Rep/Store with the following information:**
 - **Gateway IMEI # to assign the new phone number**
 - **Manufacturer: Ositech**
 - **Device Name: Titan III Wireless Cellular Gateway 4G**
- 3. Your SIM card will require certain account features activated for your Gateway to function properly. Please request the Carrier representative for the following features:**
 - **Telemetry/M2M/Machine to Machine line**
 - **Activated on at least 2-3MB data plan**
 - **Voice restriction**
 - **No static IP**
 - **Provisioned on APN vzwinternet**
- 4. To install the SIM card in your cellular Gateway, please refer to the User Manual.**

For activation on LTE networks such as AT&T or others:

- 1. Find the International Mobile Equipment Identifier (IMEI) for your cellular Gateway. For your convenience, this can be found on the peelable portion of the box label, as well as on the label applied to the device. The IMEI # is made up of 15 decimal digits.**
- 2. Contact your local AT&T Rep/Store with the following information:**
 - **Gateway IMEI # to assign the new phone number.**
 - **Manufacturer: Ositech**
 - **Model Number: T3G15401**
- 3. Your SIM card will require certain account features activated for your Gateway to function properly. Please request the Carrier representative for the following features:**
 - **Telemetry/M2M/Machine to Machine line**
 - **Activated on at least 2-3MB data plan**
 - **Voice restriction**
 - **No static IP**
 - **Provisioned on APN nxtgenphone (For AT&T Network only)**
- 4. To install the SIM card in your cellular Gateway, please refer to the User Manual.**

Once you have successfully completed the activation with your cellular Carrier, contact your LIFENET System implementation Specialist to complete the configuration and registration of your Gateway with the LIFENET System.

If there are any questions or concerns with the process of activating your Gateway(s) with your cellular Carrier, please contact LIFENET System support team at 1-800-732-3081.